

Purpose

This application note describes how an ECG (ElectroCardioGram) can be measured using a SmartPhone with a case containing Plessey Semiconductors' Electric Potential Integrated Circuit (EPIC) sensor.

Introduction

EPIC is an electrometer capable of sensing ECG signals through insulated sensors in contact with the skin. The sensors are dry-contact, so that the gels or other contact-enhancing substances normally associated with wet-electrode ECG pads are not necessary. As well as offering exciting possibilities for simplified ECG monitoring by medical professionals, this technology also makes it possible for individuals to view and collect their own detailed ECG signals on a portable device such as a SmartPhone.

SmartPhone Application

The ECG trace ideally requires two electrical signals from parts of the body on opposite sides of the heart. By mounting two sensor electrodes on the rear of a SmartPhone case, these signals are easily obtained from fingers on both hands just by holding the phone, as shown in figure 1.




Figure 1: Diagrammatic representation showing user holding smartphone and touching one sensor with each hand

To produce the ECG trace shown in figure 2 the ECG sensor in a SmartPhone application requires

1. A case - into which is built the EPIC sensors and some electronics to amplify, filter and digitise the signal and send the signal to the phone.
2. A software app to receive, process and analyse the data and display the waveform.




Figure 2: ECG trace on a smartphone

Signal processing

The collected signals should be filtered, differentially amplified and digitised by circuitry within the SmartPhone case to produce the ECG signal. Full ECG generally requires a bandwidth of 50mHz to 150Hz.

The signal is sent to the phone by bluetooth. Further analysis of the ECG trace can be performed by software with the phone app, for instance to display heart rate or other key parameters from the PQRST waveform.

Further development of the software could enable the data to be sent to – for instance – a clinic or doctor's surgery for monitoring by healthcare professionals.

For further information about this and other products, please visit:
www.plesseysemiconductors.com

Legal Notice

Product information provided by Plessey Semiconductors Limited ("Plessey") in this document is believed to be correct and accurate. Plessey reserves the right to change/correct the specifications and other data or information relating to products without notice but Plessey accepts no liability for errors that may appear in this document, howsoever occurring, or liability arising from the use or application of any information or data provided herein. Neither the supply of such information, nor the purchase or use of products conveys any licence or permission under patent, copyright, trademark or other intellectual property right of Plessey or third parties.

Products sold by Plessey are subject to its standard Terms and Conditions of Sale that are available on request. No warranty is given that products do not infringe the intellectual property rights of third parties, and furthermore, the use of products in certain ways or in combination with Plessey, or non-Plessey furnished equipments/components may infringe intellectual property rights of Plessey.

The purpose of this document is to provide information only and it may not be used, applied or reproduced (in whole or in part) for any purpose nor be taken as a representation relating to the products in question. No warranty or guarantee express or implied is made concerning the capability, performance or suitability of any product, and information concerning possible applications or methods of use is provided for guidance only and not as a recommendation. The user is solely responsible for determining the performance and suitability of the product in any application and checking that any specification or data it seeks to rely on has not been superseded.

Products are intended for normal commercial applications. For applications requiring unusual environmental requirements, extended temperature range, or high reliability capability (e.g. military, or medical applications), special processing/testing/conditions of sale may be available on application to Plessey.

Application Note 291474 issue 1

Plessey Semiconductors Ltd.

Design & Technology Centre, Delta 500, Delta Business Park, Great Western Way, Swindon, UK SN5 7XE

Tel: +44 1793 518000

Fax: +44 1793 518030

Web: www.plesseysemi.com